

## High-Stress Retention Nitrile Examination Gloves\* Accelerator Free

An alternative solution of high elasticity nitrile gloves at a lower cost

The retention feature that is commonly sought after. It is now available in nitrile at an improved performance and greater savings.

Unlike regular nitrile gloves that quickly loosen and relax after wear, our gloves maintain a skin-tight fit for an extended period.

\*Patent Pending

### KEY BENEFITS

- ✓ **Free of Sulphur & Accelerators**  
Prevent Type IV allergy
- ✓ **Highly Elastic**  
Provide a comfortable wearing experience
- ✓ **Soft & Comfortable**  
Prevent hand and finger fatigue

### PRODUCT ATTRIBUTES



Thin Gauge Design



Ambidextrous



Chlorinated



Fingertip Textured



Powder Free



Standard Cuff



Medium Tiffany Blue Colour



Not made with natural rubber latex

	ASTM Requirement			EN Requirement	
<b>Physical Dimensions</b>					
Glove Length (mm)	≥230			≥240	
<b>Thickness (mm)</b>					
Cuff	N/A			N/A	
Palm	0.07 ± 0.02			0.07 ± 0.02	
Finger	0.10 ± 0.02			0.10 ± 0.02	
Plam Width(mm) ± 4mm	Extra-Small	Small	Medium	Large	Extra-Large
	76	86	98	107	115
<b>Physical Properties</b>					
<b>Tensile (Mpa)</b>					
Before Aging	≥18			N/A	
After Aging	≥16			N/A	
<b>Elongation (%)</b>					
Before Aging	≥500			N/A	
After Aging	≥400			N/A	
<b>Force at Break (N)</b>					
Before Aging	N/A			≥6	
After Aging	N/A			≥6	

Regulatory testing is in progress and target completion by March 2022.



### CONTACT US:

<https://hartalega.com.my/contact-us/sales-contact/>

[hartalega.com.my](https://hartalega.com.my)