



























EXAMINATION GLOVES

GLOVE TYPE	WEIGHT (gram)	PALM THICKNESS		FINGER THICKNESS		COLOUR	LENGTH (mm)		SURFACE FINISH
		mil	mm	mil	mm		ASTM	EN	
Nitrile Powder Free	2.7	2.0	0.05	3.1	0.08	Blue (Dark/Ocean/Translucent) & White 	≥ 230	Median ≥ 240	Fingertips Textured
	3.0	2.2	0.05	3.1	0.08	Blue (Dark/Violet/Ocean) & White 			
	3.0	2.2	0.05	3.2	0.08	Black 			
	3.2	2.5	0.07	3.5	0.09	Blue (Violet), White & Pink 			
	3.5	3.0	0.07	3.9	0.10	Blue (Violet), White, Pink Green & Black 			
	4.0	3.2	0.08	4.3	0.11	Blue (Aqua) 			
	4.1	3.2	0.08	4.3	0.11	Orange 			
	4.5	3.4	0.08	4.7	0.12	Black 			
	4.7	3.5	0.09	4.7	0.12	Blue (Aqua) & White 			
	5.0	3.5	0.09	5.1	0.13	Blue (Aqua) 			
	5.5	4.0	0.10	5.5	0.14	Blue (Aqua) 			
	6.2	4.5	0.12	6.0	0.16	Blue (Dark) 			
Nitrile Powdered	3.5	3.0	0.07	3.9	0.10	Blue (Aqua) & White 			
Nitrile Powder Free Semi Long	3.8	2.5	0.07	4.7	0.12	White 	≥ 265		
Nitrile Powder Free Long Cuff	5.0	3.2	0.08	4.3	0.12	Blue (Violet) 	≥ 280	Median ≥ 280	
	5.8	3.5	0.09	4.7	0.12	Blue (Cobalt) 			
	6.2	4.0	0.09	5.5	0.14	Blue (Aqua) 			
	7.1	4.5	0.12	6.3	0.16	Blue (Aqua) 			
	7.8	4.8	0.13	6.7	0.17	Orange 			
Latex Powdered	5.0	4.0	0.08	4.7	0.12	White (Milky) 	≥ 230	Median ≥ 240	Smooth
	6.0	5.0	0.12	5.1	0.13	White (Milky) 			
Latex Powder Free Chlorinated	5.0	4.0	0.10	4.7	0.12	White (Milky) 	≥ 230	Median ≥ 240	Textured
	5.5	4.5	0.11	5.1	0.13	White (Milky) 			
	6.5	5.0	0.13	5.5	0.14	White (Milky) 			
Latex Polymer Coated	7.0	5.5	0.14	5.9	0.15	White (Milky) 			


SPECIALTY GLOVES
ANTIMICROBIAL GLOVES


GLOVE TYPE	WEIGHT (gram)	PALM THICKNESS		FINGER THICKNESS		COLOUR	LENGTH (mm)		SURFACE FINISH
		mil	mm	mil	mm		ASTM	EN	
Nitrile Powder Free	3.0	2.2	0.05	3.1	0.08	Blue (Violet) 	≥ 230	Median ≥ 240	Fingertips Textured
	3.2	2.5	0.07	3.5	0.09	Blue (Violet) 			
	3.5	3.0	0.07	3.9	0.10	Blue (Violet) 			

COLLOIDAL OATMEAL COATED GLOVES

GLOVE TYPE	WEIGHT (gram)	PALM THICKNESS		FINGER THICKNESS		COLOUR	LENGTH (mm)		SURFACE FINISH
		mil	mm	mil	mm		ASTM	EN	
Nitrile Powder Free COATS	3.2	2.5	0.07	3.5	0.09	Blue (Dawn) 	≥ 230	Median ≥ 240	Fingertips Textured
Latex Powder Free COATS	5.5	4.5	0.11	5.1	0.13	Green (Lemon) 			Textured


HIGH-STRESS RETENTION ACCELERATOR-FREE GLOVE

GLOVE TYPE	WEIGHT (gram)	PALM THICKNESS		FINGER THICKNESS		COLOUR	LENGTH (mm)		SURFACE FINISH
		mil	mm	mil	mm		ASTM	EN	
Nitrile Powder Free	3.5	3.0	0.07	4.0	0.10	Medium Tiffany Blue 	≥ 230	Median ≥ 240	Fingertips Textured


ACCELERATOR-FREE GLOVE

GLOVE TYPE	WEIGHT (gram)	PALM THICKNESS		FINGER THICKNESS		COLOUR	LENGTH (mm)		SURFACE FINISH
		mil	mm	mil	mm		ASTM	EN	
Nitrile Powder Free	3.2	2.5	0.07	3.5	0.09	Blue (Tiffany) 	≥ 230	Median ≥ 240	Fingertips Textured




BIODEGRADABLE GLOVE

GLOVE TYPE	WEIGHT (gram)	PALM THICKNESS		FINGER THICKNESS		COLOUR	LENGTH (mm)		SURFACE FINISH
		mil	mm	mil	mm		ASTM	EN	
Nitrile Powder Free	3.5	3.0	0.07	4.0	0.10	Blue (Violet) 	≥ 230	Median ≥ 240	Fingertips Textured


METAL DETECTABLE GLOVE

GLOVE TYPE	WEIGHT (gram)	PALM THICKNESS		FINGER THICKNESS		COLOUR	LENGTH (mm)		SURFACE FINISH
		mil	mm	mil	mm		ASTM	EN	
Nitrile Powder Free	4.7	3.2	0.08	5.1	0.13	Black 	≥ 230	Median ≥ 240	Fingertips Textured




INDUSTRIAL GLOVES

GLOVE TYPE	WEIGHT (gram)	PALM THICKNESS		FINGER THICKNESS		COLOUR	LENGTH (mm)		SURFACE FINISH
		mil	mm	mil	mm		ASTM	EN	
Nitrile Powder Free Semicon	7.0	4.5	0.11	5.9	0.15	White 	≥ 280	Median ≥ 280	Fingertips Textured
Nitrile Clean Room Class 100	7.0	4.5	0.11	5.9	0.15	White 			
Nitrile Chemical Resistant Long Cuff	11.8	8.0	0.20	11.0	0.28	Blue (Sky) 			

POLYCHLOROPRENE EXAMINATION GLOVE

GLOVE TYPE	WEIGHT (gram)	PALM THICKNESS		FINGER THICKNESS		COLOUR	LENGTH (mm)		SURFACE FINISH
		mil	mm	mil	mm		ASTM	EN	
Polychloroprene Powder Free	4.0	2.0	0.05	2.0	0.05	Turquoise & Black 	≥ 230	Median ≥ 240	Fingertips Textured

SURGICAL GLOVES

GLOVE TYPE	WEIGHT (gram)	PALM THICKNESS		FINGER THICKNESS		COLOUR	LENGTH (mm)		SURFACE FINISH
		mil	mm	mil	mm		ASTM	EN	
Latex Powdered Free	10.5	7.5	0.19	8.3	0.21	White (Cream) 	290 ± 10	Median 290 ± 10	Textured
Polyisoprene Powder Free (Underglove)	11.5	7.5	0.19	8.3	0.21	Blue 			
Polyisoprene Powder Free (Outerglove)	11.5	7.5	0.19	8.3	0.21	White (Milky) 			

All surgical products meet AQL 0.65.



CONTACT US:

<https://hartalega.com.my/contact-us/sales-contact/>